


What does the EYE-SYNC Smooth Pursuit report show?




- A** After completing the EYE-SYNC smooth pursuit assessment, the clinician receives a report with metrics measuring the individual's radial (or spatial) variability and tangential (or timing) variability, as well as the mean phase error of the eye compared to the target. The lower the score, the better dynamic visual orientation.
- B** The first plot shows the trajectory of the eye position during tracking of the circular motion of the target. Each blue dot represents an eye position captured by the camera.
- C** The second plot shows the variability of eye positions around the target. The green or red bar indicates good or poor mean phase error of the eye compared to the target respectively.
- D** The graph diagram shows an Individual's radial and tangential score compared to prior results.

The content provided is an information resource only and is not intended to be a substitute for professional medical advice, diagnosis, or treatment. Always seek the advice of your physician or other qualified health provider with any questions you may have regarding a medical condition.




What does the EYE-SYNC Saccade report show?



- A** After completing the EYE-SYNC Saccade assessment, the clinician receives a report with metrics measuring the individual's eye tracings between targets and a precision score. This is where the precision scores are located. The precision value reflects how well the individual is able to repeatedly land their eye position in the same place as they move from one target to another. The lower the score, the better precision.
- B** This is a visual illustration of where the eye lands in regards to the target and other eye positions during the assessment. This is a visual display of the precision (or lack of precision) metrics included in the report.
- C** Each eye position recorded by the camera is shown here in small blue dots. It is ideal for the lines to be as straight and close together as possible. This is where you will visualize and identify saccade like jerking recorded and included in the report.
- D** The graph diagram shows an individual's results as compared to prior performance.

The content provided is an information resource only and is not intended to be a substitute for professional medical advice, diagnosis, or treatment. Always seek the advice of your physician or other qualified health provider with any questions you may have regarding a medical condition.

What does the EYE-SYNC VOR report show?



- A** After completing the EYE-SYNC VOR assessment, the clinician receives a report with metrics measuring the individual's variability around the target position while the head is moving. This is where the precision scores are located. Precision X is the value to note after a horizontal VOR assessment is completed. Precision Y value is noted after vertical VOR assessment.
- B** This plot shows the eye position during fixed gaze on the target. Each navy blue dot represents an eye position captured by the camera. It is ideal to see the blue lines bunched along the X-axis and centered around the red dot in the middle (for a horizontal test) and along the Y-axis for a vertical test.
- C** The graph diagram shows an individual's results as compared to prior performance.

The content provided is an information resource only and is not intended to be a substitute for professional medical advice, diagnosis, or treatment. Always seek the advice of your physician or other qualified health provider with any questions you may have regarding a medical condition.